



TRi

TERUMO RADIAL INTERVENTION


EXPANDING TRANSRADIAL ACCESS FOR MORE PATIENTS

Glidesheath Slender[®]

Transradial Introducer Kit

Expanding TRI – For More Patients


Cumulative Frequency of Radial Artery Inner Diameter¹



Modified from Saito et al., Influence of the ratio between radial artery inner diameter and sheath outer diameter on radial artery flow after transradial coronary intervention.

Glidesheath Slender®

Supporting transradial intervention
every day



Expanding to Provide More Possibilities...

...with Radial Access



The proportion of patients with a mean radial diameter smaller than 2.52mm ranges from:
14% - 32% in men¹
27% - 60% in women¹



Smaller sheath footprint will help to minimize the chances of reduced radial flow, spasm, radial injury and delayed radial artery occlusion.¹



Complex PCI's are now an option: rotablator, bifurcation, left main, CTO, etc.

Features & Benefits



Glidesheath Slender®

INCREASE YOUR **RADIAL ACCESS OPTIONS**

Thin-wall technology
in a full size range
5, 6 and 7Fr

Maximize procedural outcomes

Perform diagnostic and interventional
procedures without upsizing to
a larger sheath


Enhance patient comfort

Optimal access also for small radial
arteries



Maximize procedural outcomes

Perform diagnostic and interventional procedures
without upsizing the sheath

- Unique thin-wall technology
 - 1 Fr size reduction in outer diameter
 - Maintaining larger inner-diameter equivalent



5, 6 and 7Fr Glidesheath Slender: "More with Less"



Enhance patient comfort


Reduced risk of damaging the radial artery

- M Coat™ hydrophilic coating


- Ease of insertion and removal
- Less penetration resistance than conventional sheaths¹

- TIF Tip Technology

- Total Integrated Fit, optimal tapering design for smooth penetration



Glidesheath Slender®



Perform diagnostic and interventional procedures without upsizing to a larger sheath

MAXIMIZE
PROCEDURAL OUTCOMES

Reduced risk of damaging the radial artery⁵

ENHANCE
PATIENT COMFORT

Smaller diameter sheaths enhance post-procedure hemostasis¹ and, in combination with a hydrophilic coating, reduce radial spasm and radial occlusion.^{2,3,4}

Glidesheath Slender® – Clinical results (6Fr)¹

Catheterization and Cardiovascular Interventions 00:00-00 (2013)

Original Studies

Initial Experience with the Glidesheath Slender for Transradial Coronary Angiography and Intervention: A Feasibility Study with Prospective Radial Ultrasound Follow-Up

Adel Aminian,* MD, Dariouch Dolatabadi, MD, Pascal Lefebvre, MD, Robert Zimmerman, MSc, Philippe Brunner, MD, Georges Michalakis, MD, and Jacques Lalmand, MD

- 114 patients
- All comorbid population, including ACS and complex PCI



- ✓ Procedural success rate 99,1%
- ✓ Low rate of radial artery occlusion (0,88%)
- ✓ Low rates of vascular complications and spasm
- ✓ Less invasive → shorter compression time

- Routine use of Glidesheath Slender for angiography & intervention is safe & feasible
- Potential to allow complex PCI while limiting local trauma to the artery

Glidesheath Slender® – Clinical results (5Fr)¹

Cardiovasc Interv and Ther
DOI 10.1007/s12928-015-0344-2

ORIGINAL ARTICLE

Safety and feasibility of the new 5 Fr Glidesheath Slender

Fuminobu Yoshimachi¹ • Ferdinand Kiemeneij² • Motomaru Masutani³ •
Takashi Matsukage⁴ • Akihiko Takahashi⁵ • Yuji Ikari¹



- 21 patients
- Cases done in Tergooi (NL), Dr. Kiemeneij
- TRA and TRI cases, ad-hoc PCI planned for all patients
- Right and left radial artery access
- Hemostasis with TR Band or bandage

- The new 5Fr GSS is safe and feasible for transradial angiography & intervention

- ✓ No case of radial artery occlusion
- ✓ No vascular complication or spasm
- ✓ Patent radial artery in all cases (100%)
- ✓ Low pain scores

Glidesheath Slender®


5, 6, 7Fr Glidesheath Slender

Plastic IV cannula

- 0.021" & 0.025" plastic and spring wires
- 20G & 22G Surflash needle
- Sheath length: 10cm, 16cm

Metallic entry needle

- 0.018" (5 & 6Fr 10cm sheaths), 0.021" & 0.025" spring wires
- 20G, 21G & 22G (5 & 6Fr sheaths) metallic needle
- Sheath length: 10cm, 16cm





Thank You

